




BIPHASIC DEFIBRILLATOR MONITOR CARDIOSTART

Interchangeable
reusable pads




Adult and Pediatric

5.7" colored LCD



High resolution

Easy Steps



Selection, Charge, Shock



CONTACT US

7831 NW 72ND AVENUE,
MEDLEY, FL 33166 USA

phone: +1 305 887 7552

e-mail: info@usdefib.com

BRAND TRUST AND CREDIBILITY

THERMAL PRINTER



PACEMAKER



AED MODE



SPO2




NON
INVASIVE
BLOOD
PRESSURE
NIBP



NIBP CUFF
(AVAILABLE
IN DIFFERENT
SIZES)



12VDC
AMBULANCE
INPUT FOR
MOBILE ICU



Features:

- Exponential truncated biphasic waveform
- Portable, transport handle
- 5.7" colored LCD – high resolution
- Up to 12 ECG derivations
- Thermal printer
- Ambulance input for mobile ICU

Optional parameters:

- Pulse Oximetry – (SPO2)
- Non invasive Blood Pressure (NIBP)
- Automated defibrillation with voice and text command (AED MODE)
- External Pace Maker
- Capnography (ETCO2)
- 360 Joules

